
---

# AMEDA® GLO

WEARABLE BREAST PUMP

## Instructions for Use

**Trusted** by Hospitals, **Loved** by Moms.

Ameda.com | ParentCare@ameda.com | 1.866.992.6332

---

## Table of Contents

1 Overview .....	1
2 Conventions used in this manual .....	1
3 Important safety information .....	1
3.1 Indications for use .....	1
3.2 Safety precautions .....	1
3.3 Precautions .....	2
3.4 Contraindication .....	3
3.5 Usage & operating life .....	3
4 Product description .....	3
4.1 Parts overview .....	3
4.2 Description of LED display and button function .....	3
4.3 Package content .....	4
5 Getting started .....	4
5.1 Assembly .....	4
5.2 Pumping .....	5
5.2.1 To begin pumping .....	5
5.2.2 Pumping modes .....	7
5.3 Ending your pumping session .....	8
6 Charging your breast pump .....	9
6.1 Charge the device .....	9
6.2 Battery status indicator illustration .....	10
6.3 Charging times .....	10
7 Cleaning, transfer & storage .....	10
7.1 Cleaning & maintenance .....	10
7.2 Milk storage guidelines .....	12
8 Troubleshooting .....	13
9 Technical information .....	13
9.1 Specifications .....	13
9.2 Operating conditions .....	14
9.3 Transport and storage .....	14
9.4 EMC declaration .....	15
10 Disposal information .....	17
11 Symbols definition .....	18
12 Warranty .....	20

---

## 1. Overview

Thank you for purchasing this product.

The GLO wearable breast pump features a soft silicone flange and an intuitive design to enhance comfort and convenience during use. Its easy assembly and discreet nature allow you to express milk on the go, making it a practical solution for your daily pumping routine.

## 2. Conventions Used In This Manual

### WARNING

Can lead to moderate injury.

### CAUTION

Can lead to minor injury.

### NOTE

Can lead to material damage.

### TIPS

Useful or important information that is not related to safety.

Please read this user manual carefully before using your breast pump, and save the user manual for future reference.

## 3. Important Safety Information

### 3.1 Indications for use

The Glo Wearable Breast Pump is intended to express milk from a lactating woman in order to collect milk from their breasts. The device is intended for a single user.

### 3.2 Safety Precautions

### WARNINGS

- Powered breast pumps that are designed for single users should never be rented or shared.
- Only use products or parts described in this user manual.
- Close supervision is necessary when product is used near children.

- 
- Do not modify device yourself or attempt to open or repair pump or charger.
  - Pump is not provided sterile. Boil flange, silicone diaphragm, valve and collection cup prior to first use and ensure that these parts are cooled sufficiently before being assembled.
  - Do not use harsh chemicals to clean flange, silicone diaphragm, valve and collection cup.
  - Do not wash pump or immerse in water. Wipe it with a clean, damp cloth only.
  - Do not place pump in refrigerator or freezer.
  - Do not store damp components.
  - Do not thaw breast milk in the microwave or in hot water.
  - Do not refreeze thawed breast milk.
  - Plug the charger into a visible and easily accessible electrical outlet. To disconnect from power, unplug the charger from the wall outlet.
  - Do not operate the product if the charger or charging cable is damaged by being dropped or submerged in water.
  - Do not drop or insert any object into any openings on pump.
  - Do not use pump on an airplane.
  - Do not operate where aerosol spray products are used or oxygen is being administered.
  - Pump will not work when connected to charger; always disconnect pump from charger before using.
  - Do not attempt to reset pump during use.
  - Before using the device, consult your healthcare provider to determine the appropriate duration for each pumping session.
  - Pregnant women should not use the breast pump, as it may cause contractions or premature labor.
  - If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.

### 3.3 Precautions

- Ensure flange, silicone diaphragm, valve and collection cup is clean.
- Ensure flange, silicone diaphragm, valve and collection cup are securely assembled and connected to pump.
- Before each use, check all components for cracks, chips, tears, discoloration, or wear. Discontinue use if any damage is found until parts are replaced.
- If discomfort is felt, suction can be broken by inserting a finger between the breast and the flange.
- If irritation or discomfort occurs, discontinue use and consult your healthcare provider.
- If the charger casing or wiring becomes loose, separated, or frayed, stop use of the charger immediately and contact the device manufacturer.
- Stop using the device immediately if you see smoke or detect burning from the pump unit, charger, or charging cable.

- If the package appears damaged or tampered with, do not use it and contact Ameda ParentCare immediately.

### 3.4 Contraindication


No known contraindications.

### 3.5 Usage & operating life


- Pump is reusable for a single user, and has an expected operating life of 1 year.
- Reusable parts (including flanges, silicone diaphragms and valves) should be cleaned prior to use and replacement every 3 months is recommended.

## 4. Product Description

### 4.1 Parts overview

### 4.2 LED display and button function

NO.	Component Name	Function Description
1	Battery level	Indicates the battery status
2	Expression mode	Expression mode selected on the display
3	Stimulation mode	Stimulation mode selected on the display
4	Level up	Increase suction level button

5	Level down	Decrease suction level button
6	Automatic mode	Automatic mode selected on the display
7	Massage mode	Massage mode selected on the display
8	Mode selection	Mode selection button
9	Power/Pause	Power ON/OFF and pause button
10	Level display	Indicates the suction level
11	Minutes display	Indicates how many minutes the pump has been in operation
12	Seconds display	Indicates how many seconds the pump has been in operation

### 4.3 Package Contents

Pump motor (2); Flange (2); Silicone diaphragm (2); Valve (4); Collection cup (2); Bra adjustment strap (2); USB charging cable (2); (2 each) 17, 19 & 21mm soft flange inserts; Instructions for use (1).

## 5. Getting Started

### 5.1 Assembly

To assemble the pump, follow the steps below.


#### CAUTION:

- Clean the pump as shown in these instructions.
- Wash your hands well before putting together cleaned parts.
- Assemble the pump correctly as shown in these instructions.
- To prevent leaks, wipe away any residue around the edges of parts before assembling.

(1) Install the valve in the collection cup, paying attention to orientation.


(2) Place the silicone diaphragm into the collection cup and press it down to make sure the edges seal tightly.

- 
- (3) To put the flange on the collection cup, line up the groove on the flange with the groove on the cup. Then press down on the edges with both hands to make sure it fits tight.

- (4) Put the pump motor on the collection cup. Line up the notch under the pump motor with the silicone diaphragm top. Push it down to make sure it fits snugly and seals well.


## 5.2 Pumping

### 5.2.1 To begin pumping

- (1) Press the flange tightly against your breast. Make sure there are no gaps and that your nipple is in the center. Keep the pump upright; moving it may affect performance.


- (2) Press and hold (long press) the Power/Pause  button to turn the breast pump on. The screen will light up.

(3) After the pump starts, press the Mode button to switch between modes. You can choose from four modes.


(4) Press the **M** button to switch working modes. There are four working modes (Stimulation mode, Expression mode, Massage mode and Automatic mode) for you to choose.

- The  symbol lights up indicating that the device is currently in Stimulation mode.
- The  symbol lights up indicating that the device is currently in Expression mode.
- The  symbol lights up indicating that the device is currently in Automatic mode.
- The  symbol lights up indicating that the device is currently in Massage mode.


(5) When switching to Automatic mode, the device will start from the first suction level of automatic mode by default and increase by one suction level every 15 seconds until the 12th suction level. The device will continue to operate at the 12th suction level. If you press the + or - button during this process, it will change to manual mode and stay at the corresponding suction level. At this point, you can increase/decrease the suction level by pressing the + level up button or - level down button.

(6) There are 15 suction levels in each mode. You can choose the suction that makes you comfortable by pressing the "+" level up and "-" level down button.

(7) You can pause or resume the breast pump with a short press of the  button.

(8) Press and hold (long press) the  button to turn the breast pump off.

### 5.2.2 Pumping Modes

1. This pump has level memory. It remembers the last suction level you were pumping in. The pump starts in Stimulation mode for 2 minutes every time you turn it on then it automatically moves to Expression mode unless the mode is manually changed.
2. In Automatic mode, the suction level will increase automatically. If you don't press any buttons, the suction level will go up every 15 seconds until it reaches level 12. If you press the suction level increase or decrease buttons, the automatic suction level progression will stop.
3. After the pump is paused by pressing the Power/Pause button during pumping, it will automatically shut down after 15 minutes if pumping is not restarted.
4. The breast pump will turn off automatically after being on for 30 minutes.

---

## **TIPS**


- If you regularly express more than 180mL per session on each breast, stop the expression immediately when the expressed milk reaches the 180mL level, and pour the breast milk into the bottle to prevent overfilling and overflow.
- Never tilt the breast pump when the milk has reached the maximum line on the collection cup to avoid spillage.
- The included bra adjustment straps may be attached to your bra to increase the length of your bra straps. This will give your bra enough room to hold the breast pump.

## **WARNING**


- The breast pump does not turn on while charging.

### 5.3 Ending your pumping session

- (1) When you are done expressing, press and hold (long press) the Power/Pause  button to turn the breast pump off. Break the seal between your breast and the pump flange with your finger, then carefully remove the pump.

- (2) Separate the collection cup from the pump motor.

- (3) Pour the milk into the storage bottle or bag, making sure the pour opening is lined up with the bottle or bag opening.


- (4) Clean the pump parts following the cleaning instructions.

## 6. Charging Your Breast Pump

### 6.1 Charge the device

Check the battery indicator after each use by pressing and holding (long press) the Power/Pause button. If the battery indicator blinks, charge the pump with the charger and charging cable.

1. Plug charging cable into the pump's charging port.
2. Plug other end of charging cable into the charger and plug the charger into any standard electrical outlet.
3. Once pump is charged, remove charging cable from pump and unplug charger from electrical outlet.

## 6.2 Battery status indicator illustration

status indicator displayed	Remaining pump time left
Three bars always on	~ 100 minutes
One bar - flashing	~ 5 minutes

## 6.3 Charging times

Full charge takes approximately 120 minutes and power indicator will show three bars always on when fully charged.

### TIPS

- The device has a rechargeable lithium battery inside and can only be used without the charging cable connected.
- The device has a three bar battery display. When the power is off or the charger is not connected, the battery indicator DOES NOT illuminate.
- If the charger is connected when the power is off, the battery indicator will illuminate. When the battery is charging, the battery indicator will be flashing. When the battery is fully charged, the battery indicator will be illuminated.
- When the device is operating, the battery indicator is always on when the battery is full. If the battery is low, the battery indicator will flash.
- If the charger is connected when the device is operating, the device will stop operating in order to protect the device.

### WARNING

- Use only a recommended certified 5VDC 1A adapter that meets the requirements of standard IEC 60950; and the adapter should be placed outside of the patient's environment when charging.

## 7. Cleaning, Transfer & Storage

### 7.1 Cleaning & maintenance

Pump outer housing	Wake up the pump by pressing and holding (long press) the power button.	Wipe with a clean damp cloth with 70% ethyl alcohol or isopropyl alcohol as needed. Do not immerse Pump in water.
Flange, Silicone Diaphragm, Valve, Soft flange inserts, and Collection cup	Boil in water for 10 minutes.	Clean with mild dish washing soap and clean water.

#### WARNINGS:

- Be careful when removing and cleaning the valve, if the valve is damaged, the breast pump will not work properly.
- To clean the valve, please rub it gently in the warm water with a few drops of cleaning liquid. Do not insert any objects into the valve, which may cause damage.

#### (1) Cleaning by hand

- Place pump parts in a clean wash basin used only for infant feeding items. Do not place pump parts directly in the sink.
- Add soap and hot water to basin.
- Scrub items using a clean brush used only for infant feeding items.
- Rinse by holding items under running water or by submerging in fresh water in a separate basin.
- Air-dry thoroughly. Place pump parts, wash basin, and bottle brush on a clean, unused dish towel or paper towel in an area protected from dirt and dust. Do not use a dish towel to rub or pat items dry.
- Clean wash basin and bottle brush. Rinse them well and allow them to air-dry after each use. Wash them by hand or in a dishwasher at least every few days.

#### (2) Cleaning in dishwasher

- When cleaning pump parts in a dishwasher, be sure to place small items in a closed-top basket or mesh laundry bag. Run the dishwasher using hot water and a heated drying cycle.
- Remove from dishwasher with clean hands. If items are not completely dry, place items on a clean, unused dish towel or paper towel to air-dry thoroughly before storing. Do not use a dish towel to rub or pat items dry.

#### (3) Disinfection

- Parts that come into contact with breast/breast milk, including the flange, soft flange inserts, silicone diaphragm, valve and collection cup should be disinfected by boiling in water for 10 minutes as described below.
- Separate all parts that come in contact with breast/breast milk. Wash hands thoroughly.
- Fill pot with enough water to cover all parts, bring water to a boil. Place parts in boiling water for 10 minutes.
- Allow water to cool and gently remove parts from water with tongs.
- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts in a clean, cool place when not in use. Do NOT store wet or damp parts.

## TIPS

- Flange, soft flange inserts, silicone diaphragm, valve and collection cup are dishwasher safe using top rack only.

## WARNINGS

- Do not wash pump. Wipe with clean, damp cloth, if needed.
- Use mild dish washing soap rather than harsh chemicals to clean flange, soft flange inserts, silicone diaphragm, valve, and the collection cup.
- Do not store wet or damp parts.

## 7.2 Milk storage guidelines

- (1) Write the date on the storage container. Include your child's name if you are giving the milk to a child care provider.
- (2) Gently swirl the container to mix the cream part of the breastmilk that may rise to the top back into the rest of the milk. Do not shake the milk as this can make some of the milk's valuable parts break down.
- (3) Refrigerate or chill milk right after it is pumped, if possible. You can put it in the refrigerator, place it in a cooler or insulated cooler pack, or freeze it in small (2- to 4-ounce) batches for later feedings. Pumped milk is OK without being refrigerated for up to 4 hours after pumping.

	Room Temperature 77°F (25°C)	Time in Refrigerator 39°F (4°C)	Time in Freezer 0°F (-18°C)
Freshly Expressed or Pumped	Up to 4 Hours	Up to 4 Days	6 Months - Best 12 Months - Acceptable
*Thawed, Previously Frozen	1-2 Hours	Up to 1 Day (24 Hours)	Do Not Refreeze
Leftover From Feeding (baby did not finish bottle)	Use within 2 hours after the baby is finished feeding		

\*NOTE: When thawing frozen milk, label as thawed when completely thawed (i.e. no ice crystals present). Use the time when completely thawed to base acceptable time limits for use rather than when it is taken from the freezer. These guidelines for storage and thawing of breast milk are a recommendation. Contact your lactation consultant or breastfeeding specialist for further information.

These guidelines are for healthy full-term babies and may vary for premature or sick babies. Check with your health care provider. Find more breastfeeding resources at: [WICBreastfeeding.fns.usda.gov](http://WICBreastfeeding.fns.usda.gov) and [www.cdc.gov/breastfeeding](http://www.cdc.gov/breastfeeding).

### SAFE THAWING OF BREAST MILK

Always thaw the oldest breast milk first. Remember first in, first out. Over time, the quality of breast milk can decrease.

There are several ways to thaw your breast milk:

- In the refrigerator overnight.
- Set in a container of warm or lukewarm water.
- Under lukewarm running water.
- Never thaw or heat breast milk in a microwave. Microwaving can destroy nutrients in breast milk and create hot spots, which can burn a baby's mouth.
- If you thaw breast milk in the refrigerator, use it within 24 hours. Start counting the 24 hours when the breast milk is completely thawed, not from the time when you took it out of the freezer.
- Once breast milk is brought to room temperature or warmed, use it within 2 hours.
- Never refreeze breast milk after it has thawed.

Additional information:

<https://www.cdc.gov/breastfeeding/breast-milk-preparation-and-storage/handling-breastmilk.html>

## 8. Troubleshooting

This section summarizes the most common problems you could encounter with the device. If you are unable to solve the problem with the information below, please contact the local service center in your country or contact Ameda Parent Care at 1.866.992.6332 or [parentcare@ameda.com](mailto:parentcare@ameda.com).

Problem	To solve issues
Low suction	Check that all components are assembled correctly and that the valve is positioned correctly.
	Make sure that the diaphragm, and the valve are fitted securely to the pump body and a complete seal is created.
	The air pressure in your region is too low, this may affect the suction.
No milk being expressed	Make sure the pump is assembled correctly, and the suction is created.
	Relax and try again, as expression may improve with practice.
Pain in the breast area when expressing	You may be pumping too hard. You do not need to use all the suction that the pump can generate. Only use the comfortable levels. Press the level down button to reduce the suction strength.
	Consult your health professional/breastfeeding advisor.
Leakage of breast milk during use	Remove and refit the Flange. Make sure that it is installed correctly and firmly.
	Try leaning slightly forward.

## 9. Technical Information

### 9.1 Specifications

Power requirements	
Only use the Charger recommended as below:	
Power supply	Charging voltage: 5VDC by Adapter Working voltage: 3.7VDC Li-ion battery (4W, 1200mAh each unit)
Charger specifications:	Input: 100-240 VAC, 50-60Hz, 0.4A-0.2A Output: 5VDC, 1A
Technology specifications	
Use environment	Home use
Location for use	OTC
Pump dimensions	124 mm × 111 mm × 69 mm

Package dimensions	180mm x 160mm x 140mm
Weight	252g
Collection cup capacity	180mL
Single/Multiple user device	Single user device
Provided non-sterile	Yes
Re-usable	Yes
Pump options	Single and dual
Pump type	Diaphragm pump
Stimulation vacuum range*	40 to 170mmHg
Expression vacuum range*	120 to 245mmHg
Automatic vacuum range*	40 to 245mmHg
Massage vacuum range*	40 to 140mmHg
Stimulation cycle speed*	77 to 127 cycles/minute
Expression cycle speed*	32 to 92 cycles/minute
Automatic cycle speed*	32 to 136 cycles/minute
Massage cycle speed*	79 to 143 cycles/minute
Adjustable suction levels	15 levels
Visual indicator	LED screen
Auto-off time	Automatically turns off if it is on for more than 30 minutes.
Service Life	1 year

Note: The \* symbol denotes essential performance.

### 9.2 Operating conditions

Operating Temperature	5°C to 40°C (41°F to 104°F)
Relative Humidity	15% to 90% non-condensing
Atmospheric Pressure	90 KPa to 106 KPa

### 9.3 Transport and storage conditions

Transport and storage Temperature	-20°C to 60°C (-4°F to 140°F)
Relative Humidity	15% to 90% non-condensing
Atmospheric Pressure	90 KPa to 106 KPa

---

#### 9.4 EMC Declaration

The Wearable Breast Pump may be used at home, office, or other general areas.

Special consideration should be given to use:

- **During Medical Procedures:** Do not use the Pump during any medical procedure or during imaging. For example, do not use pump during RF ablation of a skin anomaly, x-ray, MRI, or ultrasound.
- **With other equipment:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the pump and the other equipment should be observed to verify that they are operating normally.
- **In High Electrostatic Discharge Environments:** If the air is very dry, use the pump where floors are wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. If you notice shocks when you touch things, as sometimes happens in very dry winter months, be sure to touch a metal item such as a table before handling the pump.
- **In Industrial Settings:** Avoid using the pump near arc welding equipment, radar installations, and similar industrial settings.
- **Near RF Emitters:** Keep devices that use WiFi, Bluetooth, and other intentional RF emitters at least 30 cm from the pump during operation. This includes communications equipment peripherals such as antenna cables and external antennas. Otherwise, degradation of the performance of the pump could result. If you know the power output the RF emitter, the table labeled "Electromagnetic Compatibility Compliance Statement Information" can be used for more precise information about separation distances.

#### WARNINGS

- Use of accessories other than those specified by manufacturer may result in increased emissions or decreased immunity of the Wearable Breast Pump.
- This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
  - (1) This device may not cause harmful interference;
  - (2) This device must accept any interference received, including interference that may cause undesired operation.

**NOTE**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules.

**Electromagnetic emission**

<b>Guidance and manufacturer's declaration – electromagnetic emission</b>		
The model W8S is intended for use in the electromagnetic environment specified below. The customer or the user of the model W8S should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The model W8S uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The W8S is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

**Electromagnetic immunity**

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The model W8S are intended for use in the electromagnetic environment specified below. The customer or the user of the model W8S should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge IEC 61000-4-2	$\pm 2$ kV, $\pm 4$ kV, $\pm 6$ kV, $\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	$\pm 2$ kV, $\pm 4$ kV, $\pm 6$ kV, $\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/ burst IEC 61000-4-4	$\pm 2$ kV for power supply lines 100 kHz repetition frequency $\pm 1$ kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV differential mode line-line	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply	0 % $U_T$ (100 % dip in $U_T$ ) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % $U_T$	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the models W8S requires continued operation during power mains interruptions, it is recommended that the models W8S be powered from an uninterruptible power supply or a

input lines IEC 61000-4-11	(100 % dip in $U_T$ ) for 1 cycle at 0°  70 % $U_T$ (30 % dip in $U_T$ ) for 25/30 cycles at 0°  0 % $U_T$ (100 % dip in $U_T$ ) for 250/300 cycle at 0°		battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
IMMUNITY to proximity magnetic fields	65A/m, Modulation: Pulse modulation, 2.1KHz Test frequency:134.2KH Z; 7.5A/m, Modulation: Pulse modulation, 50KHz Test frequency:13.56M Hz;	65A/m, Modulation: Pulse modulation, 2.1KHz Test frequency:134.2KH Z; 7.5A/m, Modulation: Pulse modulation, 50KHz Test frequency:13.56M Hz;	/
NOTE: $U_T$ is the a. c. mains voltage prior to application of the test level.			
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands and amateur radio bands <sup>a</sup>	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands and amateur radio bands <sup>a</sup>	Portable and mobile RF communications equipment should be used no closer to any part of the models W8S, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance</b> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80MHz to 800MHz  $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800MHz to 2.7GHz where P is the maximum output

Radiated RF IEC 61000-4-3	10 V/m 80MHz to 2.7 GHz	10 V/m 80 MHz to 2.7 GHz	<p>power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m).<sup>b</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>c</sup> should be less than the compliance level in each frequency range<sup>d</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
------------------------------	----------------------------	-----------------------------	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a /The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b /The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c/ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the / is used exceeds the applicable RF compliance level above, the / should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the /.

d/ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distances

#### Recommended separation distances between portable and mobile RF communications equipment and the model W8S

The model W8S is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model W8S can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model W8S as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output of transmitter	Separation distance according to frequency of transmitter /m
-------------------------------------	--

/W	150 kHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.035	0.07
0.1	0.38	0.11	0.22
1	1.2	0.35	0.7
10	3.8	1.11	2.21
100	12	3.5	7

For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Recommended separation distances between RF wireless communications equipment

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

Frequency /MHz	Maximum Power /W	Distance /m	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
385	1.8	0.3	27	27	RF wireless communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. <b>Recommended separation distance</b>  Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance
450	2	0.3	28	28	
710	0.2	0.3	9	9	
745					
780					
810	2	0.3	28	28	
870					
930					

1720	2	0.3	28	28	in meters (m). Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 	
1845						
1970						
2450	2	0.3	28	28		
5240	0.2	0.3	9	9		
5500						
5785						
Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.						

## 10. Disposal Information

The Breast Pump contains a rechargeable lithium ion battery. Use your local return and collection system for electrical waste and electronic equipment (including batteries).

Flange, silicone diaphragm, soft flange inserts, and valve should be disposed of in your home's waste management container.

### WARNINGS

- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

## 11. Symbols Definition

The warning symbol identifies all instructions that are important for safety. Failure to observe these instructions can lead to injury or damage to the breast pump. When used in conjunction with the following words, the warning symbol stands for:

### (1) Pump control buttons


Button symbol	Function description
	ON/OFF or Pause/Start
	Level down = decrease suction level
	Level up = increase suction level
	Mode selection

### (2) Symbols on the packaging

Symbol	Definition
	This symbol indicates the manufacturer.
	This symbol indicates to keep the device dry.
	This symbol indicates to handle the fragile device with care.
	This symbol indicates to keep the device away from sunlight.
	This symbol indicates that waste electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately.
	This symbol indicates the temperature limitation for operation, transport and storage.
	This symbol indicates the humidity limitation for operation, transport and storage.
	This symbol indicates the atmospheric pressure limitation for operation, transport and storage.
	This symbol indicates that a product has been assessed by the manufacturer and deemed to meet EU safety, health and environmental protection requirements.

	This symbol indicates to follow instructions for use.
	This symbol indicates that the packaging contains products intended to come in contact with food according to regulation 1935/2004.
	This symbol indicates that the packaging is recyclable.
	This symbol indicates a carton package.

### (3) SYMBOLS ON THE DEVICE

	This symbol indicates to follow instructions for use.
	This symbol indicates the manufacturer.
	This symbol indicates that waste electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately.
	This symbol indicates a type BF applied part.
<b>IP22</b>	This symbol indicates the degree of ingress protection (IP).
	This symbol indicates the date of manufacture (four digits for the year and two digits for the month).
	This symbol indicates manufacturer's serial number of the device.
	This symbol indicates direct current.
	This symbol indicates that a product has been assessed by the manufacturer and deemed to meet EU safety, health and environmental protection requirements.
	This symbol indicates that the device's electromagnetic interference is below the limit set by the US Federal Communications Commission.

This device conforms to the following standards:

AAMI/IEC 60601-1 Medical Electrical Equipment - General requirements for basic safety and essential performance

IEC 60601-1-2 Medical Electrical Equipment - Electromagnetic disturbances

IEC 60601-1-11 Medical Electrical Equipment - Home healthcare environment

The manufacturer is ISO 13485 certified.

## 12. Warranty

**THE AMEDA WARRANTY IS NONTRANSFERABLE BETWEEN INDIVIDUALS. THE AMEDA GLO ELECTRIC BREAST PUMP ("PRODUCT") IS WARRANTED TO THE ORIGINAL USER – WHETHER INSTITUTIONAL OR INDIVIDUAL ("FIRST USER") – ONLY.**

From the date of the Individual User's purchase of this product, Ameda, Inc. ("Ameda") warrants the Product to the Individual User against defects in material or workmanship for a period of one year for the pump mechanism and 90 days on the included partsexcept for breast pump diaphragms and valves. Ameda's sole obligation under this limited warranty shall be to replace, at Ameda's option, any Product that is determined to be defective by Ameda and determined to be covered by this express limited warranty. Replacement under this limited warranty is the sole andexclusive remedy of the First User.

**PROOF OF PURCHASE IN THE FORM OF A RECEIPTED INVOICE OR BILL OF SALE EVIDENCING THAT THE PRODUCT IS WITHIN THE WARRANTY PERIOD MUST BE PRESENTED TO OBTAIN WARRANTY SERVICE.** This limited warranty is extended by Ameda ONLY to the First User and is not assignable or transferable. For warranty assistance contact Ameda at **1.866.992.6332**.

EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ON THE PRODUCT IS HEREBY DISCLAIMED. AMEDA SHALL NOT BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES OF ANY KIND FOR BREACH OF THE EXPRESS LIMITED WARRANTY ON THE PRODUCT OR ANY WARRANTY IMPLIED BY OPERATION OF LAW. OTHER THAN THE EXPRESS LIMITED WARRANTY SET FORTH ABOVE, THERE ARE NO OTHER WARRANTIES THAT ACCOMPANY THIS PRODUCT AND ANY ORAL, WRITTEN OR ANY OTHER REPRESENTATION OF ANY KIND IS HEREBY DISCLAIMED.

### NON-APPLICABILITY OF WARRANTY

In addition to any other limitations on the warranty discussed above, the warranty shall expire and be of no force or effect upon the occurrence of any of the following:

- (i) the cause of damage to the Product is due to improper handling or use.  
opening of the pump motor, modifications, self-repair or use outside the instructions in this manual voids warranty.
- (iii) the Product is not properly maintained by the Individual User as determined by Ameda.
- (iiii) non-Ameda approved accessories or spare parts are used with the Product as determined by Ameda.

If any of the above occurs, Ameda shall not be responsible for damage to the Product, parts thereof or injuries arising therefrom, either directly or indirectly. The warranty terms set forth above replaces all prior warranties with respect to the purchased Product, whether in writing or otherwise.

**Outside the United States:** Call your local distributor or location where you purchased the Product. For a listing of distributors in your country, please visit [www.ameda.com](http://www.ameda.com).

 Manufactured by: Shenzhen TPH Technology Co., Ltd.

Distributed worldwide by:  
Ameda, Inc.  
8206 Cross Park Drive  
Austin, TX 78754  
United States  
1.866.992.6332  
[www.ameda.com](http://www.ameda.com)

Distributed in Canada by:  
Distribué au Canada par:  
Mothers Choice Products Ltd.  
2133-11871 Horseshoe Way  
Richmond, BC V7A5H5  
1.800.604.6225  
[www.motherschoiceproducts.com](http://www.motherschoiceproducts.com)

Distributed in Australia by:  
Distribué en Australie par:  
Distribuido en Australia por:  
Midmed Pty Ltd.  
Unit 4, 62 Borthwick Avenue  
61 (7) 3348 9155  
[www.midmed.com.au](http://www.midmed.com.au)


Pump 1 Serial No:																					
Pump 2 Serial No:																					

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