
AMEDA® GLO

WEARABLE BREAST PUMP

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Trusted by Hospitals, **Loved** by Moms.

Ameda.com | ParentCare@ameda.com | 1.866.992.6332

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Electromagnetic emission

| Guidance and manufacturer's declaration – electromagnetic emission | | |
|---|-------------------|--|
| The model W8S is intended for use in the electromagnetic environment specified below. The customer or the user of the model W8S should assure that it is used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment – guidance |
| RF emissions CISPR 11 | Group 1 | The model W8S uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The W8S is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions CISPR 11 | Class B | |
| Harmonic emissions IEC 61000-3-2 | Not applicable | |
| Voltage fluctuations / flicker emissions IEC 61000-3-3 | Not applicable | |

Electromagnetic immunity

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|--|--|---|---|
| The model W8S are intended for use in the electromagnetic environment specified below. The customer or the user of the model W8S should assure that it is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment -guidance |
| Electrostatic discharge IEC 61000-4-2 | ±2 kV, ±4 kV, ±6 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air | ±2 kV, ±4 kV, ±6 kV, ±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrostatic transient/ burst IEC 61000-4-4 | ± 2 kV for power supply lines 100 kHz repetition frequency ± 1 kV for input/output lines | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 0.5 kV, ± 1 kV differential mode line-line | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply | 0 % U_T (100 % dip in U_T) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U_T | Not applicable | Mains power quality should be that of a typical commercial or hospital environment. If the user of the models W8S requires continued operation during power mains interruptions, it is recommended that the models W8S be powered from an uninterruptible power supply or a |

| | | | |
|--|--|---|---|
| input lines IEC 61000-4-11 | (100 % dip in U_T) for 1 cycle at 0° 70 % U_T (30 % dip in U_T) for 25/30 cycles at 0° 0 % U_T (100 % dip in U_T) for 250/300 cycle at 0° | | battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m, 50/60Hz | 30 A/m, 50/60Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| IMMUNITY to proximity magnetic fields | 65A/m, Modulation: Pulse modulation, 2.1KHz Test frequency:134.2KH z; 7.5A/m, Modulation: Pulse modulation, 50KHz Test frequency:13.56M Hz; | 65A/m, Modulation: Pulse modulation, 2.1KHz Test frequency:134.2KH z; 7.5A/m, Modulation: Pulse modulation, 50KHz Test frequency:13.56M Hz; | / |
| NOTE: U_T is the a. c. mains voltage prior to application of the test level. | | | |
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands and amateur radio bands ^a | 3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands and amateur radio bands ^a | Portable and mobile RF communications equipment should be used no closer to any part of the models W8S, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.7\text{GHz}$ where P is the maximum output |

| | | | |
|------------------------------|----------------------------|-----------------------------|--|
| Radiated RF IEC 61000-4-3 | 10 V/m 80MHz to 2.7 GHz | 10 V/m 80 MHz to 2.7 GHz | <p>power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^c should be less than the compliance level in each frequency range^d</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
|------------------------------|----------------------------|-----------------------------|--|

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a /The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b /The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

c/ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the / is used exceeds the applicable RF compliance level above, the / should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the /.

d/ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances


| Recommended separation distances between portable and mobile RF communications equipment and the model W8S | |
|---|--|
| The model W8S is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model W8S can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model W8S as recommended below, according to the maximum output power of the communications equipment. | |
| Rated maximum output of transmitter | Separation distance according to frequency of transmitter /m |

| /W | 150 kHz to 80 MHz $d = [\frac{3.5}{V_1}] \sqrt{P}$ | 80 MHz to 800 MHz $d = [\frac{3.5}{E_1}] \sqrt{P}$ | 800 MHz to 2.7 GHz $d = [\frac{7}{E_1}] \sqrt{P}$ |
|------|---|---|--|
| 0.01 | 0.12 | 0.035 | 0.07 |
| 0.1 | 0.38 | 0.11 | 0.22 |
| 1 | 1.2 | 0.35 | 0.7 |
| 10 | 3.8 | 1.11 | 2.21 |
| 100 | 12 | 3.5 | 7 |

For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

| Recommended separation distances between RF wireless communications equipment | | | | | |
|--|------------------|-------------|----------------------|------------------|---|
| The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment. | | | | | |
| Frequency /MHz | Maximum Power /W | Distance /m | IEC 60601 Test Level | Compliance Level | Electromagnetic Environment – Guidance |
| 385 | 1.8 | 0.3 | 27 | 27 | RF wireless communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance  |
| 450 | 2 | 0.3 | 28 | 28 | |
| 710 | 0.2 | 0.3 | 9 | 9 | |
| 745 | | | | | |
| 780 | | | | | |
| 810 | 2 | 0.3 | 28 | 28 | |
| 870 | | | | | |
| 930 | | | | | |

| | | | | | |
|--|-----|-----|----|----|---|
| 1720 | 2 | 0.3 | 28 | 28 | <p>in meters (m). Field strengths from fixed RF transmitter, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| 1845 | | | | | |
| 1970 | | | | | |
| 2450 | 2 | 0.3 | 28 | 28 | |
| 5240 | 0.2 | 0.3 | 9 | 9 | |
| 5500 | | | | | |
| 5785 | | | | | |
| <p>Note 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | | | |

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Valmistaja: Shenzhen TPH Technology Co., Ltd.
Markkinoija: Ameda Inc., 8206 Cross Park Drive, Austin, TX 78754, United States, 1.866.992.6332